

# Koala *Adult Bacterial Viral Filter*



**Code: GM-KO-39016**

**Features:**

1. Equipped with 0.3 micron bacterial filter membrane for up to 99.99% of efficiency.
2. Resistance less than 1.6cmH<sub>2</sub>O/30lpm.
3. Standard coaxial 22F/15M for inlet port and 22M/15F for patient port meet ISO 5356-1 standard, 4.3 luer connector with strap for pressure sensing or CO<sub>2</sub> detection.
4. Material: Polypropylene.



**Code: GM-KO-39013**

**Features:**

1. Equipped with 0.3 micron bacterial filter membrane for up to 99.99% of efficiency.
2. Pressure resistance less than 1.6cmH<sub>2</sub>O/30lpm.
3. Standard coaxial 22F/15M for inlet port and 22M/15F for patient port meet ISO 5356-1 standard.
4. Material: Polypropylene.

Distributed by:



**Koala Medical Pty Ltd**

PHONE: 02 9979 9374

FAX: 02 9979 6409

Email: [KOALAMEDICAL@bigpond.com](mailto:KOALAMEDICAL@bigpond.com)

[www.koalamedical.com.au](http://www.koalamedical.com.au)