

Koala *Autoclavable Swivel Connectors and Extension Sets*



The products are made of high grade plastics and silicones formulated to withstand the rigors of cleaning, processing and autoclaving without affecting their performance.

SWIVEL CONNECTORS

- Clear for visualization of secretions, obstructions and humidity
- 360° Swivel connections with ISO tapers and O-ring seals

SUCTION/BRONCHOSCOPY PORT


- Permits access without disruption of ventilation
- Functions as either a cap seal or plug seal
- Sealaround port allows for PEEP maintenance during procedure

EXTENSION SETS

- Smooth Inner Diameter
- Low Dead Space
- Minimal Resistance
- Flexes with patient/circuit repositioning
- Rigid machine and connector for security

100% LATEX FREE

DHD60-0010 – Each
Double Swivel Elbow with Suction Port



DHD60-1510 – 5/box
Silicone Extension Tube with 15mm inlet and double swivel elbow with suction port

DHD60-0002 – Each
Single Swivel Elbow without Port



DHD60-2210 – 5/box
Silicone Extension Tube with 22mm inlet and double swivel elbow with suction port

DHD60-0003 – Each
Straight Swivel 15mm OD to 22mm OD fixed/15mm ID Swivel



DHD60-1502 – 5/box
Silicone Extension Tube with 15mm inlet and single swivel elbow without Port

DHD60-1503 – 5/box
Silicone Extension Tube with 15mm inlet and straight swivel (15mm OD – 22mm OD Fixed/15mm ID Swivel)

Distributed by:



Koala Medical Pty Ltd

PHONE: 02 9979 9374

FAX: 02 9979 6409

Email: KOALAMEDICAL@bigpond.com

www.koalamedical.com.au