

Koala Disposable Breathing Circuits

Drappable Tubing - Fully Welded



Adult 72" Breathing Circuit

Code

Complete 72" standard tubing with swivel wye, Blue Vent Limb 2 x 22mm Male Male connectors, 2 litre Latex Free Rebreathing Bag
Compatible with "Ulco" Anaesthetic machines, no extra adapters required.

GM-KO-6693
20 Circuits per box



Adapters for the "Datex-Ohmeda/GE Healthcare S/5 ADU" Anaesthetic machines.

Reusable Stainless Steel
22mm.

KO.S5.22MM
Each

ADU Circuit: Complete 72" standard tubing with swivel wye, Blue Vent Limb with built-in ADU adapter, 22mm Male Male connectors, 2 litre Latex Free Rebreathing Bag, 60" White four limb for connecting from the ventilator to the absorber

GM-KO-6693/60/ADU
15 Circuits per box

Distributed by:



Koala Medical Pty Ltd

PHONE: 02 9979 9374

FAX: 02 9979 6409

Email: KOALAMEDICAL@bigpond.com

www.koalamedical.com.au