

Koala *Disposable Connectors*

Connectors and Adaptors

Features

- All connections conform to ISO standard for taper and fit
- Single use,
- Latex Free
- 50 per box



GM-KO-72211
Mask Elbow
15mm ID/22mm OD
To 15mm OD



GM-KO-76112
Gas Sampling Elbow
15mm ID/22mm OD
To 15mm O.D.



GM-KO-74055
Multi Adapter
15mm ID/22mm OD
To 22mm OD



GM-KO-73046
Washington Tee
15mm ID/15mm OD



GM-KO-75322
Gas Sampling Tee
15mm ID/22mm OD
To 15mm OD



OL-001515
Gas Sampling Tee
15mm ID/22mm OD
To 15mm OD



GM-KO-73017
Adult Swivel Wye
15mm ID/22mm OD
To two 22mm OD



GM-KO-73018
**Adult Swivel Wye with
Gas Sampling Port**
15mm ID/22mm OD
To two 22mm OD



GM-KO-72235
Metered Dose Elbow
15mm ID/15mm OD

100% LATEX FREE

Distributed by:



Koala Medical Pty Ltd

PHONE: 02 9979 9374

FAX: 02 9979 6409

Email: KOALAMEDICAL@bigpond.com

www.koalamedical.com.au