

Koala Medical

Disposable Silicone Laryngeal Airway



The Koala “Colour Coded” Disposable Silicone Sterile Laryngeal Mask Airway and Reinforced Laryngeal Mask Airway is Highly Comparable to the Original Laryngeal Mask Airway with Epiglottis Bars.

A Complete Range of Colour Coded Sizes are Available in both “classic” and reinforced in Neonate to Large Adult Sizes.

Disposable Silicone Laryngeal Airway

DESCRIPTION	SIZE	COLOUR	WEIGHT (KG)	CUFF FILL VOLUME (ML)	CATALOGUE
NeoNatal	1	Black	> 5 kg	4	LMA1004S2
Paediatric	1.5	White	5-10 Kg	7	LMA1507S2
Infant	2	Red	10-20 Kg	10	LMA2010S2
Child	2.5	Orange	20-30 Kg	14	LMA2514S2
Large Child/ Small Adult	3	Green	30-50 Kg	20	LMA3020S2
Adult	4	Yellow	50-70 Kg	30	LMA4030S2
Large Adult	5	Purple	70 Kg>	40	LMA5040S2

Disposable Silicone Reinforced Laryngeal Airway

DESCRIPTION	SIZE	COLOUR	WEIGHT (KG)	CUFF FILL VOLUME (ML)	CATALOGUE
Neonatal	1	Black	>5 Kg	4	LMARF1004S2
Paediatric	1.5	White	5-10 Kg	7	LMARF1504S2
Infant	2	Red	10-20Kg	10	LMARF204S2
Child	2.5	Orange	20-30 Kg	14	LMARF2514S2
Large Child/ Small Adult	3	Green	30-50 Kg	20	LMARF3020S2
Adult	4	Yellow	50-70 Kg	30	LMARF4030S2
Large Adult	5	Purple	70-100 Kg	40	LMARF5040S2



KOALA
MEDICAL PTY LTD

Phone:

02 9979 9374

Fax:

02 9979 6409

Email:

KOALAMEDICAL@bigpond.com