


# Koala Disposable Breathing Circuits

## Expandable Tubing - Fully Welded



**The unique Koala Extend-Flex Circuit is available in both Adult and Paediatric. Circuit expands, bends and flexes to desired length and direction. NO GAS LEAKS**

### Adult 22mm x 110"

#### Code

2 x 22mm expandable hoses to approximately 110", 1 x Parallel-Wye, 1 x Mask Elbow without Luer Port, 1 x 22mm Blue Vent Limb expandable to 110", 1 x 2.0 Litre Rebreathing Bag

**GM-KO-67017+2**  
**20 Circuits per box**

### Paediatric 15mm x 110"

#### Code

2 x 15mm expandable hoses to approximately 110", 1 x Parallel-Wye, 1 x Mask Elbow without Luer Port

**GM-KO-67016**  
**40 Circuits per box**

2 x 15mm expandable hoses to approximately 110", 1 x Parallel-Wye, 1 x Mask Elbow without Luer Port with 1 x 15mm Blue Vent Limb expandable to 110" and 1 x Rebreathing Bag 0.5 Litre

**GM-KO-67016+15MM-0.5**  
**20 Circuits per box**

2 x 15mm expandable hoses to approximately 110", 1 x Parallel-Wye, 1 x Mask Elbow without Luer Port with 1 x 15mm Blue Vent Limb expandable to 110" and 1 x Rebreathing Bag 1 Litre

**GM-KO-67016+15MM-1.0**  
**20 Circuits per box**

Distributed by:



**Koala Medical Pty Ltd**

PHONE: 02 9979 9374

FAX: 02 9979 6409

Email: [KOALAMEDICAL@bigpond.com](mailto:KOALAMEDICAL@bigpond.com)

[www.koalamedical.com.au](http://www.koalamedical.com.au)