

"Koala Blue"

100% Neoprene Anaesthetic Face Mask



Information

The Colour Coded Neoprene Masks are unique due to their design and choice of material. The use of Neoprene means the masks are latex-free and therefore ideal for patients and staff who are allergic to latex. Neoprene is less reactive to external factors, making the life-span of these masks longer than that of comparable products.

Profiled Inflatable Cuff to facilitate optimum air-seal with patient's face makes this an ideal anaesthetic mask to replace (old-fashioned) Black Rubber Masks.

Decontamination /Disinfection by Washer/Disinfector as per Australian Standards AS4187

Complete with one-way valve for easy cuff inflation/deflation

	Size 1 – Infant	Code # DI131000
	Size 2 – Child Small	Code # DI132000
	Size 3 – Child Medium / Large	Code # DI133000
	Size 4 – Adult Small	Code # DI134000
	Size 5 – Adult Medium	Code # DI135000
	Size 6 – Adult Large	Code # DI136000
	One-way Valves	Code # DIV2300V

Distributed by:



Koala Medical Pty Ltd

PHONE: 02 9979 9374

FAX: 02 9979 6409

Email: KOALAMEDICAL@bigpond.com

www.koalamedical.com.au