

Koala Resuscitators

Disposable



- PVC Manual resuscitator
- PC Reservoir valve set
- PC intake valve set
- Air cushion face mask
- Oxygen tubing 2 metres
- PE Reservoir bag
- POP-OFF Valve
- LATEX FREE

Reusable



- Silicone Manual resuscitator
- Autoclavable up to 134°C
- PSF Reservoir valve set
- PSF intake valve set
- Silicone Air cushion face mask
- Oxygen tubing 2 metres
- PVC Reservoir bag
- POP-OFF Valve
- LATEX FREE

<u>Size & Description</u>	<u>Code Disposable</u>	<u>Code Reusable</u>
Adult - 1500ml bag, 2500ml Reservoir bag, Size 5 Mask, 60cmH ₂ O pop-off	60101	60301
Child - 550ml bag, 2500ml Reservoir bag, Size 3 Mask, 40cmH ₂ O pop-off	60102	60302
Infant - 280ml bag, 600ml Reservoir bag, Size 1 Mask, 40cmH ₂ O pop-off	60103	60303

Distributed by:



Koala Medical Pty Ltd

PHONE: 02 9979 9374

FAX: 02 9979 6409

Email: KOALAMEDICAL@bigpond.com

www.koalamedical.com.au