


# Koala Reusable Silicone *Laryngeal Airway*



## *Information*

The Koala “*Cost Effective*” Laryngeal Airway is comparable to the original laryngeal mask airway without the epiglottis bars. The product therefore can provide a conduit for tracheal intubation, especially in a difficult airway situation. Clear access for suction catheter or fiberscope through the un-obstructed lumen is possible.

The Silicone Laryngeal Airway is recommended for 40 uses, with individual tracking provided with each airway.

A complete range of sizes are available from neonate to large adult, with smaller sizes in half size increments.

<b>Description</b>	<b>Size</b>	<b>Recommended Weight (kg)</b>	<b>Cuff Fill Volume (ml)</b>	<b>Code</b>
Neonatal	1.0	> 5 kg	4	2-A101-1.0
Paediatric	1.5	5-10 kg	7	2-A102-1.5
Infant	2.0	10-20 kg	10	2-A103-2.0
Child	2.5	20-30 kg	14	2-A104-2.5
Large Child/Small Adult	3.0	30-50 kg	20	2-A105-3.0
Adult	4.0	50-70 kg	30	2-A106-4.0
Large Adult	5.0	70 > kg	40	2-A107-5.0

*Multiple purchases attracts discount*

Distributed by:



**Koala Medical Pty Ltd**

PHONE: 02 9979 9374

FAX: 02 9979 6409

Email: [KOALAMEDICAL@bigpond.com](mailto:KOALAMEDICAL@bigpond.com)

[www.koalamedical.com.au](http://www.koalamedical.com.au)