

Super LED Laryngoscope Handle


Provides a light intensity source many lumens greater than available devices with technology well beyond existing laryngoscopes handles.

“Seeing it, is believing it”



Super LED Laryngoscope Medium Handle

Code : SG-01-2600



Premium Quality

- Super Light Emitting Diode.
- Super Intensity High Performing Light Source.
- 10000 Hour Diode Life.
- Quality All Brass Microswitch.
- Brushed Stainless Steel Finish.
- New Ergonomic Shape Ensures Comfort and Grip.
- ISO compatible with All Green Spec Fibre Optic Laryngoscope Blades.



- All Brass Microswitch with Light Source



KOALA
MEDICAL PTY LTD

Phone: 02 9979 9374
Fax: 02 9979 6409
Email: KOALAMEDICAL@bigpond.com