



Disposable Adult Anaesthetic Circuits for Low Flow Anaesthesia



Adult Anaesthetic Circuits

- Complete Adult 72" or 120" standard tubing with swivel wye
- Blue vent limb
- 2 X 22mm Male/Male connectors
- Free 2.0 Litre rebreathing bag
- Compatible with "Ulco" Anaesthetic Machines
- No extra adaptors required

KO.S5.22MM

- Adaptors for the "Datex Ohmeda" and "GE Healthcare S/5 ADU" Anaesthetic Machines
- Reusable stainless steel 22mm

ITEM	CODE	ORDER QUANTITY
Adult 72" Circuit	GM-KO-6693	Box of 20
Adult 120" Circuit	GM-KO-6693/120	Box of 10
Adaptors for "Datex Ohmeda" & GE Healthcare machines	KO.S5.22MM	Each

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au







Disposable Paediatric Circuits

Ayres T-Piece with Rotary Scavenging Valve

Code: GM-KO-6400

- Blue scavenging tube
- 1.5m fresh gas line
- 0.5 Litre rebreathing bag
- Box of 15



Rotary Scavenging Valve

Ayres T-Piece Circuit without Scavenging

Code: GM-KO-6401

- 1.5M fresh gas line
- 0.5 Litre rebreathing bag
- Box of 15



DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au



Joey Circuit



- 15 mm or 10 mm diameter breathing tubes
- 15 mm or 10mm diameter blue ventilation limb
- 2 X 22mm Male/Male connectors
- 10mm Circuit supplied with 0.5L rebreathing bag
- 15mm Circuit supplied with 1.0L rebreathing bag
- Compatible with “ULCO” Anaesthetic Machines with no extra adaptors required



ITEM	CODE
Medisize Hygrovent Child Angled	300 510 000
Medisize Hygrovent Child Straight	300 520 000

ITEM	DESCRIPTION	CODE	ORDER QUANTITY
Paediatric Circuit 72” 15mm	Complete 72” 15mm Paediatric Circuit with 1.0 Litre re-breathing bag	GM-KO-67000	Box of 20
Paediatric Circuit 120” 15mm	Complete 120” 15mm Paediatric Circuit with 1.0 Litre re-breathing bag	GM-KO-67000/120	Box of 15
Neonatal Circuit 72” 10mm	Complete 72” 10mm Neonatal Circuit with 0.5 Litre re-breathing bag	GM-KO-67010	Box of 20
Adaptor	Adaptor for “Datex-Ohmeda” & “GE Healthcare S/5 ADU”	KO.S5.22M	Each

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au





Septum Co-Axial Circuit

Centre Divided Single Limb Circuit



- Single limb circuit capable of low flow anaesthesia
- 3m length
- Light weight design aids in airway management
- Heat exchanged between limbs provides warmer inspired gas

ITEM	DESCRIPTION	CODE	ORDER QUANTITY
Septum Circuit 1.8M	Coaxial Septum Circuit 1.8m length with blue extendable rebreathing limb and 2.0 Litre rebreathing bag	GM-KO-46658	Box of 20
Septum Circuit 3.0M	Coaxial Septum Circuit 3.0m length with blue extendable rebreathing limb and 2.0 Litre rebreathing bag	GM-KO-46657	Box of 20



Ultrabain – Disposable Coaxial System, with APL Valve

The Armstrong Medical Ayre's Tee Piece Circuit incorporates the ULTRA APL Valve as the means of pressure control and release within the system



- Suitable for spontaneous and controlled breathing
- Designed to increase flexibility around the patient's head and neck area
- Ideal for the induction and maintenance of anaesthesia during ophthalmic, dental and laryngeal surgery
- Lightweight tubing minimises drag on patient's airway
- Scavenging adapter can be removed if necessary

APL VALVE Features

- Maximum pressure limit – 60cmH₂O
- Low flow resistance at open position – less than 1cmH₂O
- Unique scavenging system
- Audible cap rotation with intimate movement
- Scavenged gases prevented from re-entering breathing system
- Performance not affected by moisture & humidity in respiratory gases

ITEM	CODE	ORDER QUANTITY
Ultrabain Coaxial Circuit	AMCX4606/6/1LF	Box of 20

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au





Disposable Breathing Circuits Expandable Tubing – Fully Welded



The unique Koala Extend-Flex Circuit is available in both Adult and Paediatric. Circuit expands, bends and flexes to desired length and direction. NO GAS LEAKS

ITEM	DESCRIPTION	CODE	ORDER QUANTITY
Adult 22mm x 110"	2 x 22mm expandable hoses to approx 110", 1 x Parallel-Wye, 1 x Mask Elbow without Luer Port, 1 x 22mm Blue Vent Limb expandable to approx 110", 1 x 2.0 Litre Rebreathing Bag	GM-KO-67017+2	Box of 20
Paediatric 15mm x 110"	2 x 15mm expandable hoses to approx 110", 1 x Parallel-Wye, 1 x Mask Elbow without Luer Port	GM-KO-67016	Box of 20
Paediatric 15mm x 110" (0.5 Litre)	2 x 15mm expandable hoses to approx 110", 1 x Parallel-Wye, 1 x Mask Elbow without Luer Port, 1 x 15mm Blue Vent Limb expandable to approx 110" and 1 x Rebreathing Bag (0.5 Litre)	GM-KO- 67016+15MM-0.5 (0.5 Litre)	Box of 20
Paediatric 15mm x 110" (1.0 Litre)	2 x 15mm expandable hoses to approx 110", 1 x Parallel-Wye, 1 x Mask Elbow without Luer Port, 1 x 15mm Blue Vent Limb expandable to 110" and 1 x Rebreathing Bag (1.0 Litre)	GM-KO- 67016+15MM-1.0 (1.0 Litre)	Box of 20



KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
*Specialising in Patient Monitoring • Laryngoscopes
Airway Management • Anaesthetic Devices
Infant Care • Maternity Care and more...*

**100%
LATEX
FREE**

Smooth Bor Tubing & Paediatric Circle Circuit



ITEM	SIZE	CODE
Smooth Bor Adult Tubing 22mm Female/Female Cuffs	4" (10cm)	SB403-4
	12" (30cm)	SB403-12
	24" (61cm)	SB403-24
	48" (122cm)	GM403-48
	60" (152cm)	GM403-60
	72" (183cm)	GM403-72
	96" (236cm)	GM403-96
Smooth Bor Paediatric Tubing 15mm Female/Female	12" (30cm)	SB405-12FF
	24" (61cm)	SB405-24FF
	60" (152cm)	SB405-60FF

ITEM	SIZE	CODE
Smooth Bor Paediatric Tubing 15mm Male/Male	12" (30cm)	SB405-12MM
	24" (61cm)	SB405-24MM
Smooth Bor Paediatric Tubing 15mm Male/Female	12" (30cm)	SB405-12MF
	24" (61cm)	SB405-42MF
	60" (152cm)	SB405-60MF
Smooth Bor Vent Tubing 22mm Female/ Female BLUE Cuffs	60" (152cm)	GM403-60-BC
	72" (183cm)	GM403-72-BC
	96" (236cm)	SB403-96-BC

ITEM	DESCRIPTION	CODE
Paediatric Circle Circuit	Smooth Bor Tubing 60"FF x 2 Smooth Bor Tubing 60"MF x 1 Silicone Rebreathing Bag 0.5 Litre Paediatric Fixed Wye Machine Cuffs x 3 Straight Bag Adapter x 1	KOA405103-60

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au

