




**100%
LATEX
FREE**

Reusable Connectors



REUSABLE TUBE INHALER CODE = DHD60-0139

SWIVEL CONNECTOR FOR THE SPRAYING OF
AEROSOL MEDICATION INTO TRACHEAL TUBES.


- No Interruption in Respiration.
- With Suction Port.
- Swivel Connection.
- Easy to Use.
- Auto Clavable.



DHD60-0002

ITEM (A) = Autoclavable, (R) = Reusable	CODE
Adapter (A) 22mm M/15mm F to 22mm M/15mm F	GM74043
Adapter (A) 15mm M/10mm F to 15mm M/15mm F	900MR533
Adapter (A) 22mm F/22mm M to 22mm M/15mm F	GM74044
Adapter (A) 15mm M/10mm F to 22mm F/22mm M	GM74042
Adapter (A) 15mm M/10mm F to 22mm M/15mm F	GM74041

ITEM (A) = Autoclavable, (R) = Reusable	CODE
Single Swivel Elbow Only (A) 15mm ID/22mm OD-15mm OD	DHD60-0002
Elbow Only (R) 15mm ID/22mm OD-15mm OD	OL-018901
Elbow with Luer Port (R) 15mm ID/22mm OD-15mm OD	OL-018901L



ITEM (A) = Autoclavable, (R) = Reusable	CODE
Paediatric Fixed Wye (A) 22mm M/15mm F to 15mm M	GM73037
Adult Fixed Wye (A) 22mm M/15mm F to 22mm M	GM73001
Adult Swivel Wye (A) 22mm M/15mm F to 22mm M	DHD60-1400

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au



Single Use Connectors



Code: GM-KO-73017 & GM-KO-73018

- Single use swivel wye feature 2 x 22mm
- Limb connection and 22/15 mm patient
- Connector ISO fittings and are available with and without Luer Port



Code: GM-KO-72015 and GM-KO-72111

- Single use mask elbow feature 15mm
- Male and 15mm Female connections
- ISO fittings and are available with and without Luer Port



Code: GM-KO-76228

- Single use T-Piece with Luer Port 15mm
- Male and 22mm/15mm Male ISO fittings



Code: GM-KO-1423

- Single use metered dose elbow
- 22mm/15mm Male ISO fittings

ITEM	CODE	ORDER QUANTITY
Single Use Swivel Wye with Luer Port	GM-KO-73018	Box of 50
Single Use Swivel Wye without Luer Port	GM-KO-73017	Box of 50
Mask Elbow With Luer Port	GM-KO-72015	Box of 50
Mask Elbow Without Luer Port	GM-KO-72111	Box of 50
Tee-Piece with Luer Port	GM-KO-76228	Box of 50
Metered Dose Elbow	GM-KO-1423	Box of 50

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au



KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
*Specialising in Patient Monitoring • Laryngoscopes
Airway Management • Anaesthetic Devices
Infant Care • Maternity Care and more...*

**100%
LATEX
FREE**

Disposable Licorice Sticks



GM-KO-11



GM-KO-1516FF-2



GM-KO-12



GM-KO-4475FF-2

- Maximum Flexibility and Low Resistance to Air Flow
- Angled Fixation and Dead-space Control with Flexi Tubes
- All Connectors Conform with ISO Standard for Taper and Fit
- Single Patient Use, 50 per case

ITEM	DESCRIPTION	CODE
Flexi Tube Straight Hard Connectors	Expandable 15 mm x 15 cm Tube 22 mm Male/15mm Male Hard Connector	GM-KO-12
Flexi Tube Straight	Expandable 15mm x 15cm Tube 22 mm Female Hard/15mm Female Soft connector	GM-KO-11
Flexi Tube with Mask Elbow	Expandable 15 mm x 15 cm Tube 22 mm Male/15 mm Female Elbow 15 mm Male Circuit Connection	GM-KO-1518FF-2
Flexi Tube with Mask Elbow and Gas Sampling Port	Expandable 15 mm x 15 cm Tube 22 mm Male/15 mm Female Elbow 15mm Male Circuit Connection	GM-KO-1516FF-2
Flexi Tube with Double Swivel Elbow and Suction Port	Expandable 15 mm x 15 cm Tube 22 mm Male/15 mm Female Elbow 15 mm Male Hard Circuit Connection	GM-KO-4475FF-2

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au





KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
*Specialising in Patient Monitoring • Laryngoscopes
Airway Management • Anaesthetic Devices
Infant Care • Maternity Care and more...*

**100%
LATEX
FREE**

Autoclavable Swivel Connectors and Extension Sets



The products are made of high grade plastics and silicones formulated to withstand the rigors of cleaning, processing and autoclaving without affecting their performance.

SWIVEL CONNECTORS

- Clear for visualization of secretions, obstructions and humidity
- 360° Swivel connections with ISO tapers and O-ring seals

SUCTION/BRONCHOSCOPY PORT

- Permits access without disruption of ventilation
- Functions as either a cap seal or plug seal
- Sealaround port allows for PEEP maintenance during procedure

EXTENSION SETS

- Smooth Inner Diameter
- Low Dead Space
- Minimal Resistance
- Flexes with patient/circuit repositioning
- Rigid machine and connector for security

ITEM	CODE
Double Swivel Elbow with Suction Port	DHD60-0010
Single Swivel Elbow without Port	DHD60-0002
Straight Swivel 15mm OD to 22mm OD fixed/15mm ID Swivel	DHD60-0003
Silicone Extension Tube with 15mm inlet and double swivel elbow with suction port (5/box)	DHD60-1510
Silicone Extension Tube with 22mm inlet and double swivel elbow with suction port (5/box)	DHD60-2210
Silicone Extension Tube with 15mm inlet and single swivel elbow without Port (5/box)	DHD60-1502
Silicone Extension Tube with 15mm inlet and straight swivel (15mm OD – 22mm OD Fixed/15mm ID Swivel) (5/box)	DHD60-1503

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au

