



Super LED Ultra-Light Green Handle



Available in Medium, Small and Stubby Sizes

The new light weight autoclavable polymer Super LED Handle

- Slim light weight design.
- High Quality LED components.
- Autoclavable Polymer Shell.

CODE	DESCRIPTION	ORDER QUANTITY
SG-01-2600-LWM	Medium Ultra-Light Handle	Each
SG-01-2600-LWD	Stubby Ultra-Light Handle	Each
SG-01-2600-LWS	Small Ultra-Light Handle	Each

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au



KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
*Specialising in Patient Monitoring • Laryngoscopes
Airway Management • Anaesthetic Devices
Infant Care • Maternity Care and more...*

SUPER LED Laryngoscope Handle

Provides a light intensity source many lumens greater than available devices with technology well beyond existing laryngoscopes handles.

***"Seeing it is
believing it!"***


Premium Quality

- Super Light Emitting Diode.
- Super Intensity High Performing Light Source.
- 10000 Hour Diode Life.
- Quality All Brass Microswitch.
- Brushed Stainless Steel Finish.
- New Ergonomic Shape Ensures Comfort and Grip.
- ISO compatible with All Green Spec Fibre Optic Laryngoscope Blades.




Super LED Laryngoscope
Medium Handle

Code: SG-01-2600



Super LED Laryngoscope
Stubby Handle

Code: SG-01-2600-D



Super LED Laryngoscope
Small Handle

Code: SG-01-2600-S



DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409


Email: sales@koalamedical.com.au



KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
*Specialising in Patient Monitoring • Laryngoscopes
Airway Management • Anaesthetic Devices
Infant Care • Maternity Care and more...*

All Metal Single Use Fibre Optic and Conventional Laryngoscope Blades



Quality Disposable Full Metal Blades at Affordable Prices.

- Eliminate Cross-Infection
- Clean Ready to Use
- ISO Compatible with both 'Green & Conventional' Systems



ITEM	CODE
No. 0 Neonate Macintosh	SG-LL7110
No. 1 Baby Macintosh	SG-LL7111
No. 2 Child Macintosh	SG-LL7112
No. 3 Adult Macintosh	SG-LL7113
No. 4 Large Adult Macintosh	SG-LL7114
MILLER FIBRE OPTIC 10 PER BOX	CODE
No. 0 Neonate Miller	SG-LL7220
No. 1 Baby Miller	SG-LL7221
MACINTOSH CONVENTIONAL 10 PER BOX	
No. 0 Neonate Macintosh	SG-LL3110
No. 1 Baby Macintosh	SG-LL3111
No. 2 Child Macintosh	SG-LL3112
No. 3 Adult Macintosh	SG-LL3113
No. 4 Large Adult Macintosh	SG-LL3114



DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au



KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
*Specialising in Patient Monitoring • Laryngoscopes
Airway Management • Anaesthetic Devices
Infant Care • Maternity Care and more...*

Maxi-Lite Macintosh Enclosed Fibre Optic Laryngoscope Blades



Experience White Light...

the new generation Maxi-Lite range of enclosed Fibre Optic Laryngoscopes

ITEM	CODE
Macintosh No. 0 Neonate Blade	SG-LE7110
Macintosh No. 1 Baby Blade	SG-LE7111
Macintosh No. 2 Child Blade	SG-LE7112
Macintosh No. 3 Adult Blade	SG-LE7113
Macintosh No. 4 Large Adult Blade	SG-LE7114
Macintosh No. 5 Extra Large Blade	SG-LE7115
Super Led Handle	SG-01-2600

**SPECIAL
OFFER**

Purchase 5 or More
Max-Lite Blades and Receive
One Free Super-LED Handle



DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au






Fibre Optic Difficult Intubation Blades & Handles



Kessell Blade & Stubby Handle

The increased angle of the Kessell laryngoscope blade of about 110 degrees may allow easier insertion of the laryngoscope blade in obstetric patients.

ITEM	CODE
Kessell Blade Size 3	SG-KS7113H
Kessell Blade Size 4	SG-KS7114H
Stubby Handle	SG-01-2600-D



McCoy Type Laryngoscope Blade for Difficult Intubation

Koala McCoy Laryngoscopes feature an articulated distal tip to lift the epiglottis allowing greater visualization during difficult intubation.

SIZE	CODE
Size 2 Macintosh	SG-MC7112
Size 2 Macintosh with Handle	SG-MC7112-1
Size 3 Macintosh	SG-MC7113
Size 3 Macintosh with Handle	SG-MC7113-1
Size 4 Macintosh	SG-MC7114
Size 4 Macintosh with Handle	SG-MC7114-1

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au



KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
*Specialising in Patient Monitoring • Laryngoscopes
Airway Management • Anaesthetic Devices
Infant Care • Maternity Care and more...*

Maxi-Lite Miller Enclosed Fibre Optic Laryngoscope Blades



Experience White Light... the new generation Maxi-Lite range of Fibre Optic Laryngoscopes with enclosed fibre guide.

ITEM	CODE
Miller No. 00 Premmie Blade	SG-LE7229
Miller No. 0 Neonate Blade	SG-LE7220
Miller No. 1 Baby Blade	SG-LE7221
Miller No. 2 Child Blade	SG-LE7222
Miller No. 3 Adult Blade	SG-LE7223



DISTRIBUTED BY

Koala Medical Pty Ltd


Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au



Video Laryngoscope with Reusable Blades



HM-RMAC Blades



HM-VL3R

- Ergonomic lightweight design, anti-microbial handle.
- High resolution video function with internal 4G SD Memory Card, able to take photos and record videos for training/education
- 3.5 inch full colour, non-glare display, convenient for users to observe the laryngeal structure
- 2M pixel excellent Optics Camera, anti-fog design with six LED light, CMOS camera chip
- Micro USB port – recharging or file uploading
- 316 medical stainless steel reusable blade, easy for sterilization by hydrogen peroxide gas plasma or cleaner, 1000 times
- Various size blades available

ITEM	CODE	ORDER QUANTITY
Reusable Video Laryngoscope with 3 Reusable Blades (MAC 2, MAC 3, MAC 4)	HM-VL3R	Each
Reusable Blade MAC 1	HM-RMAC1	Each
Reusable Blade MAC 2	HM-RMAC2	Each
Reusable Blade MAC 3	HM-RMAC3	Each
Reusable Blade MAC 4	HM-RMAC4	Each
Reusable Blade MAC D5	HM-RMACD5	Each
Reusable Blade MILLER 00	HM-RMIL00	Each
Reusable Blade MILLER 0	HM-RMIL0	Each
Reusable Blade MILLER 1	HM-RMIL1	Each

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au





KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
Specialising in Anaesthetic Devices &
Airway Management – Laproscopic
Surgery – Maternity & Infant Care

Video Laryngoscope with Disposable Blades



HM-VL3 DII



HM-DMAC Blades

- Ergonomic lightweight design, anti-microbial handle.
- High resolution video function with internal 4G SD Memory Card, able to take photos and record videos for training/education
- 3.5 inch full colour, non-glare display, convenient for users to observe the laryngeal structure
- 2M pixel excellent Optics Camera, anti-fog design with six LED light, CMOS camera chip
- Micro USB port – recharging or file uploading
- Various size blades available

ITEM	CODE	ORDER QUANTITY
Disposable Video Laryngoscope with 3 Disposable Blades (MAC 2, MAC 3, MAC 4)	HM-VL3D DII	Each
Disposable Blade MAC 1	HM-DMAC1 DII	Pack of 5
Disposable Blade MAC 2	HM-DMAC2 DII	Pack of 5
Disposable Blade MAC 3	HM-DMAC3 DII	Pack of 5
Disposable Blade MAC 4	HM-DMAC4 DII	Pack of 5
Disposable Blade MAC 5	HM-DMAC5 DII	Pack of 5



DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au