



KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
*Specialising in Patient Monitoring • Laryngoscopes
Airway Management • Anaesthetic Devices
Infant Care • Maternity Care and more...*

**BRAND
NEW**

Enhanced Nitrous Circuits

Now Longer with the Added Benefit of a Directional Safety Connector that Directs Gas Away From the Midwives Face



GM-KO-6727 now expandable to 450cm

For every 1 circuit purchased KOALA MEDICAL will supply 1 Umbilical Cord Clamp (Sterile) AT NO EXTRA COST eg: BOX of 30 supplied with 30 Cord Clamps

**SPECIAL
OFFER**

Purchase 600* Circuits
and receive a FREE Doppler

*600 Obstetric Circuits must be purchased over a 12 month period or as one bulk order. **Conditions Apply**.

ITEM	DESCRIPTION	CODE
Blue Obstetric Circuit 30 per Box	Expandable Blue Flexi Tube from 122cm to 450cm. Tee Piece with One Way Valve. Mouth Piece. Bacterial Viral Filter.	GM-KO-6727
Scavenging Obstetric Circuit 30 per Box	Expandable Blue Flexi Tube from 122cm (48") to 3.4 m (134"). Pink Scavenging Flexi Tube from 122cm (48") to 3.4 m (134"). Tee Piece with One Way Valve. Mouth Piece. Bacterial Viral Filter.	GM-KO-6727SCAV

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au





Sterile Umbilical Cord Clamps and Cord Clamp Cutters



Clamp Cutters

- Sterile
- Single use
- Quick and easy to use
- Light and ergonomic

Cord Clamps

- Sterile
- Single use
- Yellow in colour

CODE	DESCRIPTION	ORDER QUANTITY
PAH-CUTTER	Umbilical Cord Clamp Cutters	Box of 50
PAH-003007	Umbilical Cord Clamps – Single Clamp	Pack of 100

**SPECIAL
OFFER**

Cord Clamps **FREE** of Charge when
Obstetric Circuits are purchased



KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
Specialising in Anaesthetic Devices &
Airway Management – Laproscopic
Surgery – Maternity & Infant Care

Baby Eye Shades



Y-B Type Eye Shade



Y-C Type Eyeshade



Y-B Type Eye Shade

ITEM	CODE	ORDER QUANTITY
Protective Eye Shades Y-B Type – Small	OT-YB-S	Pack of 20
Protective Eye Shades Y-B Type – Medium	OT-YB-M	Pack of 20
Protective Eye Shades Y-B Type – Large	OT-YB-L	Pack of 20
Protective Eye Shades Y-C Type – Small	OT-YC-S	Pack of 20
Protective Eye Shades Y-C Type – Medium	OT-YC-M	Pack of 20
Protective Eye Shades Y-C Type – Large	OT-YC-L	Pack of 20

- Made from high quality nonwoven fabric elastic material with special contour design
- Soft, comfortable and safe
- Suitable for various baby's head shapes, available in two designs
- Available in three sizes – Small, Medium & Large

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au





Baby Arm Boards



- Newborn transfusion hand support
- Specially designed, pre-shaped board provides maximum immobilization on the infant's limb when securing IV lines
- Unique "wings" combined with the bendable yet firm internal aluminum support provides an ideal fit
- Soft fabric foam cushion improves comfort and protection
- Baby boards are latex free
- Available in three sizes – Small, Medium and Large

ITEM	CODE	ORDER QUANTITY
Arm Boards – Small 1.5 x 6.5cm	OT-AB-S	Pack of 20
Arm Boards – Medium 2.5 x 9.8cm	OT-AB-M	Pack of 20
Arm Boards – Large 3 x 12cm	OT-AB-L	Pack of 20



DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au



KOALA
Medical Pty Ltd

WORLD CLASS MEDICAL SUPPLIES,
Specialising in Anaesthetic Devices &
Airway Management – Laproscopic
Surgery – Maternity & Infant Care

Novos Bililed Maxi+ Phototherapy Unit



Novos Bililed Maxi+ is a new advanced phototherapy unit with 50,000 hours of LED performance & extended effective surface area. Novos Bililed Maxi+ high intensity phototherapy treats infants with hyperbilirubinemia.



- Optimum wavelength. Blue Super LED's emit light focused on between 440-460 nm region so bilirubin breakdown is maximized.
- 16 pcs super LEDs provide uniform light distribution over the infant's whole body surface area.
- Provides an extended effective surface area to cover whole body surface of even active term infants.
- Light intensity is adjustable at 5 different levels according to each patient's need.
- Removable and lightweight head unit to use on incubators especially in constricted NICU's.
- Head can be angled to use with warmers.
- Adjustable height to use together with baby cribs, incubators or warmers.
- Easy maneuverability of castors to slide easily under warmers, incubators and cribs.
- Technical specifications available.



CODE	DESCRIPTION	BOX QUANTITY
NOV-037000	Bililed Maxi+ Phototherapy Unit	1
NOV-037013	Bililed Maxi+ Mobile Stand	1

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au

