



**KOALA**  
Medical Pty Ltd

**WORLD CLASS MEDICAL SUPPLIES,**  
*Specialising in Patient Monitoring • Laryngoscopes  
Airway Management • Anaesthetic Devices  
Infant Care • Maternity Care and more...*

## GiO Endotracheal Tube Pressure Gauge

Accurate, precise and easy to use



ITEM	COLOUR	CODE
0-3cm H <sub>2</sub> O	PINK	GM-K-GB10
0-80cm H <sub>2</sub> O	BLUE	GM-K-GB20
-30-100cm H <sub>2</sub> O	PURPLE	GM-K-GB30
-60-0cm H <sub>2</sub> O	GREEN	GM-K-GB40

The Koala GiO Digital Pressure Gauge allows for fine monitoring in an ICU & operating theatre environment. It is simple and easy to adjust air pressure in the ETT cuff with a high degree of accuracy and precision. Adjustments can be by 0.1 cmH<sub>2</sub>O in a measurement range is from 0 to 80cmH<sub>2</sub>O.

Monitoring can be done whenever needed, both during and after procedures.

- Backlight to enable reading in a dark environment.
- Real time display of pressure.
- Displays peak pressure reading when P&H function is activated.
- No need for pump interaction to check cmH<sub>2</sub>O pressure.
- Minimal intervention.
- Simple, compact and easy to use.

### DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)





**KOALA**  
Medical Pty Ltd

**WORLD CLASS MEDICAL SUPPLIES,**  
*Specialising in Patient Monitoring • Laryngoscopes  
Airway Management • Anaesthetic Devices  
Infant Care • Maternity Care and more...*

## H100B Pulse Oximeter



Neonate Wrap Probe Available Separately.

ITEM	CODE
Pulse Oximeter with Silicone Boot and Adult Probe	ED H100B
Adult Pulse Oximeter Probe	ED-H100B Adult Probe
Neonate Wrap Pulse Oximeter Probe	ED-H100B Neonate Wrap
Pole Mounting Kit	ED-H100B Mount Kit

Compact Pulse Oximeter ideal for spot checks and transportation. Comes Complete with Adult Probe and Silicone Protective Boot.

The H100B is also Pole Mountable when using the Mount Kit.

Neonatal Mode with Separate Wrap Style Probe Available.

- LCD Display
- Pulse Tone
- Audible and Visible Alarm
- Long Battery Life up-to 50 hours using 4 AA batteries
- Backlight Control Automatic Power Off for Power Saving
- SpO<sub>2</sub> and Pulse Rate Measurement
- Display Numeric Value and Pleth Waveform Simultaneously
- Data Storage Capacity which can be downloaded to a PC
- 10 Minute Review of SpO<sub>2</sub> and Pulse Rate Display Function
- Compatible with Nellcor (OXIMAX) and BCI DB9 sensors



### DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)





# Therm-O-Finder



## NON-CONTACT IR THERMOMETER

- Reading from Temporal Artery
- Quick Accurate Non-Contact Temperature Detection
- No covers required with No Risk of Virus Transmission
- Back-light Display for Night Use
- Auto Power Off

## SPECIFICATION

- Accuracy  $\pm 0.2^{\circ}\text{C}$
- Operation Temperature  $15-40^{\circ}\text{C}$
- Measurement Distance 2-3cm from forehead center
- Measurement Range Body Temperature  $34.0-42.5^{\circ}\text{C}$
- Response Time = 1 second



ITEM	CODE
ThermOfinder Non Contact Thermometer	FS-300

## DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)