



## Sigmoidoscope



**The Evexar Self-Illuminating Sigmoidoscope** is a single use, fully disposable device that utilises a convenient (universal) Illuminator inserted into the base of the handle, and a light rod to direct light to the end of the scope where it is most needed.

This enables examination and therapeutic procedures to be performed in clinics and the operating room, giving complete clinical freedom of use.

Sigmoidoscopes and Illuminators are completely disposable as a unit, after use.

- Ideal for Specialist Clinicians for use in Hospitals and Private Clinics, as well as for General Practitioners in their surgeries
- Self-illuminated; no need for fibre optic cables or cleaning
- Integrated obturator design prevents obstruction of the examination channel and light guide
- Fully disposable as a unit after use
- Reduced risk of cross-infection since the scope is single use and disposable
- Convenient and cost effective

CODE	DESCRIPTION	ORDER QUANTITY
EV-01111	Sigmoidoscope with Illuminator	Box of 25

### DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)







## Proctoscope 19mm and 25mm

**The Evexar Self-Illuminating 19mm and 25mm Proctoscopes** are single use, fully disposable devices that utilise a convenient (universal) illuminator inserted into the base of the handle, and a light rod to direct light to the end of the scope where it is most needed.

This enables examinations and minor therapeutic procedures to be performed in all areas such as Out-Patient Departments, Accident and Emergency, General Practitioner surgeries, as well as in the operating room, giving complete clinical freedom to use anywhere.



Proctoscopes and Illuminators are completely disposable as a unit, after use.

- Ideal for General Practitioner Surgeries and Out - Patient Clinics; wherever examinations are performed
- Ideal for Specialist Clinicians for use in Hospitals and Private Clinics
- Self-illuminated; no need for fibre optic cables or cleaning
- Integrated obturator design prevents obstruction of the examination channel and light guide
- Fully disposable as a unit after use
- Reduced risk of cross-infection since the scope is single use and disposable
- Convenient and cost effective

CODE	DESCRIPTION	ORDER QUANTITY
EV-01107	19mm Proctoscope with Illuminator	Box of 25
EV-01112	25mm Proctoscope with Illuminator	Box of 25

**DISTRIBUTED BY**

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)





# Anchor Grip with Gel Fixation and SoftFix™ Balloon



## Features:

- Visible Trocar Design with universal seal
- Accommodates a wide range of instruments
- Detachable cannula cap for easy tissue removal
- Allows easy single hand adjustment of cannula depth and the ability to suture the cannula in place
- Obturator with click in design for secure and stable trocar placement
- Provides atraumatic and air-tight external fixation
- Softfix™ latex free balloon enhances internal stabilization and minimizes fascial trauma

CODE	DESCRIPTION	ORDER QUANTITY
UN-FBT105100	Balloon Trocar 5mm/10cm	Box of 5
UN-FBT112100	Balloon Trocar 12mm/10cm	Box of 5

## DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)







# Tropian Tech Filtering Trocar



## Features:

- Triple Layer Microporous Filter
- Filter purifies smoke and removes hazardous material from air
- Micro Pores 1 - 50µm can absorb aeroform toxic material
- Small in size however filter has large surface area of 1200 m<sup>2</sup> / 1g
- Dual seal airshields with minimal CO<sub>2</sub> leakage during procedure
- Bladeless Tip Trocar
- Transparent cannula allows enhanced visualisation



CO<sub>2</sub> leak escaping from top of filter is purified gas



CODE	DESCRIPTION	ORDER QUANTITY
TT-ZT-22F	Trocar with Triple Layer Filter 5mm/10cm	Box of 10

## DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)







## Tropian Tech Standard Trocar



### Features:

- Good grip ergonomic designed Trocar palm
- Ribbed cannula provides stable holding position
- ZT-30 and ZT-40 integrated converter allowing the use of 5-12mm diameter instruments
- Dual seal airshields with minimal CO2 leakage during procedure
- Bladeless Tip Trocar
- Transparent cannula allows enhanced visualisation
- Detachable sample collection port



### TT-ZT-22F

Trocar with Triple Layer Filter 5mm/10cm  
(refer separate brochure)

CODE	DESCRIPTION	ORDER QUANTITY
TT-ZT-10	Trocar 5mm/7cm	Box of 10
TT-ZT-20	Trocar 5mm/10cm	Box of 10
TT-ZT-30	Trocar 11mm/10cm	Box of 10
TT-ZT-40	Trocar 12mm/10cm	Box of 10

### DISTRIBUTED BY

**Koala Medical Pty Ltd**


Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)



# Unimax Visible Surgical Trocars



Clear cannula for increased visualisation. Clear tip allows direct sight of abdominal wall layers

## Unimax Visible Trocar:

- Provides safe entry into the peritoneal cavity
- Reduces trocar related injuries
- Reduces scarring
- No fascial closure is required
- Less pain and herniation rate
- Ribbed cannula for optimal attachment
- Ergonomic handle allows easy insertion
- Detachable clear cannula head enables easier specimen retrieval
- Universal seal to accommodate 5 - 12 mm instruments

## Auto Locking Bladed Trocar's also Available

- Locking Button for Safety Shield Trocar Mechanism.
- Available in Kit Form  
1 Trocar 2 Ports



CODE	DESCRIPTION	ORDER QUANTITY
UN-FVT705070	Visible Trocar 5mm/7cm	Box of 5
UN-FVT705100	Visible Trocar 5mm/10cm	Box of 5
UN-FVT710100	Visible Trocar 10mm/10cm	Box of 5
UN-FVT712100	Visible Trocar 12mm/10cm	Box of 5
UN-FLK300300	Visible Trocar 1 x 5mm/10cm Trocar 2 x 5mm/10cm Cannulas 1 x 10mm/10cm Hasson Trocar	Box of 5
UN-FLK301300	Visible Trocar 1 x 5mm/10cm Trocar 2 x 5mm/10cm Cannulas 1 x 10mm/10cm Visible Trocar	Box of 5

CODE	DESCRIPTION	ORDER QUANTITY
UN-FSTD05070	5mm/7cm Trocar	Box of 5
UN-FSTD05100	5mm/10cm Trocar	Box of 5
UN-FSTD10100	10mm/10cm Trocar	Box of 5
UN-FSKD05100	1 x 5mm/10cm Trocar 2 x 5mm Cannulas	Box of 5

\* Custom Trocar Kits are also available on request

## DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)





# Surgical Light Handle Cover Universal



This disposable surgical light handle cover is universal and unique and fits most of the O.R. surgical lights

## Features:

- Convenient and cost-effective
- Quick and easy operation
- Reduce set up time and the cost of sterilization of the handle
- Pre-compressed state of package saves shipping cost and storage space



CODE	DESCRIPTION	ORDER QUANTITY
PAH-09086	Universal Light Handle Cover	Box of 50

## DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)





# Laparoscopic Monopolar Scissors



## Laparoscopic Monopolar Scissors:

- Stainless steel blades sustain sharpness throughout procedures
- Precise and sharp blades for smooth cutting
- Extended PTFE insulation reduces the chance of thermal damage
- 360 degree wheel allows the shaft to rotate 360 degrees
- Ergonomic handle increases comfort during operation and provides effortless operation

CODE	DESCRIPTION	ORDER QUANTITY
UN-FMS533050	Monopolar Scissors 5mm/33cm (long tip)	Box of 5
UN-FMS533054	Monopolar Micro Scissors 5mm/33cm (1/2") (short tip)	Box of 5

## DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409


Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)



# Laparoscopic Grasping Instruments



Complete Range of Grasping and Dissecting Instruments Available



CODE	DESCRIPTION	ORDER QUANTITY
UN-FEG533050	Straight Grasping Forcep 5mm/33cm	Box of 10
UN-FED533050	Maryland Dissecting Forcep 5mm/33cm	Box of 10
UN-FEC533050	Clinch Grasping Forcep 5mm/33cm	Box of 10
UN-FEB533050	Babcock Grasping Forcep 5mm/33cm	Box of 10
UN-FRT533050	Rat Tooth Grasping Forcep 5mm/33cm	Box of 10
UN-FDG533050	Duck Bill Grasping Forcep 5mm/33cm	Box of 10
UN-FEJ533050	Johan Grasping Forcep 5mm/33cm (38mm Long Clip)	Box of 10
UN-FEJ533051	Johan Grasping Forcep 5mm/33cm (22mm Short Clip)	Box of 10

## DISTRIBUTED BY



**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)



## PROPLATE



The advanced gel of PROPLATE will provide you with benefits of improved conductivity, lower skin resistance, residue-free and minimized risk of pad site burn.

### Optimal Adhesive Strength

With advanced technology, we improve the adhesion between skin and itself.

### Lower Impedance

With optimal Hydro-Gel, it can maintain high current density.


### Residue-Free

Well made Hydro-Gel Plate provides you with the best performance without the gel mess.

Backing Line

Gel Line

Aluminium Line



### ESU Plate / Safety Ring

#### Minimizes Risk of pad site burns

PRO PLATE ESU pads are designed and manufactured using an advanced gel overlay system to reduce the possibility of edge burn or bite

POPULATION	MODEL		TYPE	SIZE	CONDUCTIVE AREA	TOTAL AREA	ORDER QUANTITY
	Uncorded Cloth Type / Foam Type	Corded Cloth Type / Foam Type					
Adult	BP-9641C / 9641F	P9641C / P9641F	Single	154*124	115	169	Box of 100
Adult	BP-9642C / 9642F	P9642C / P9642F	Split (Dual)	154*124	108	169	Box of 100
Pediatric	BP-9631C / 9631F	P9631C / P9631F	Single	135*124	99	148	Box of 100
Pediatric	BP-9632C / 9632F	P9632C / P9632F	Split (Dual)	135*124	92	148	Box of 100

#### DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)








# PROPENCIL

The PROPENCIL, Electrosurgical Unit Pencil, is available with button switch for smooth activation.

## TELESCOPIC SMOKE PENCIL

### The easiest way to extend the length

- Still be able to focus on the surgery area whilst the length is extended
- Taking only 1 second to have the adequate length
- The easiest way to pull and push back the telescopic suction tube with blade
- Thin and light design gives the best comfortability to the users
- Best suction capability
- Extremely flexible air hose



### Release Button

The telescope can be varied from 2.5 to 6.0 inches with just a push of the release button and sliding the electrode tube front and back.

MODEL	BUTTON TYPE	PRODUCT DESCRIPTION	CABLE LENGTH	ORDER QUANTITY
BP-EvBN	Push Button	2.5" Stainless Coated Blade, No Holster	10" Cable	Box of 10
BP-EvBNH	Push Button	2.5" Stainless Coated Blade, Holster	10" Cable	Box of 10

### DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)





## Needle Holders

Pistol  $\Phi 5 \times 330\text{mm}$



"V" Shape  $\Phi 5 \times 330\text{mm}$



CE 0197



Straight



Curved



Self-righting



Self-righting



Self-righting

- Medical Grade Stainless Steel
- Strong, reliable grip
- Ergonomic, smooth control
- "A" Shape Needle Holders available on request

ITEM	CODE	ORDER QUANTITY
Needle Holder, "Pistol" Straight, $\Phi 5 \times 330\text{mm}$	HK-101.102	Each
Needle Holder, "Pistol" Curved, $\Phi 5 \times 330\text{mm}$	HK-101.102A	Each
Needle Holder, "Pistol" Self-Righting, $\Phi 5 \times 330\text{mm}$	HK-101.102B	Each
Needle Holder, "Pistol" Self-Righting, $\Phi 5 \times 330\text{mm}$	HK-101.102C	Each
Needle Holder, "Pistol" Self-Righting, $\Phi 5 \times 330\text{mm}$	HK-101.102D	Each
Needle Holder, "V" Shape Straight, $\Phi 5 \times 330\text{mm}$	HK-101.024	Each
Needle Holder, "V" Shape Curved, $\Phi 5 \times 330\text{mm}$	HK-101.024A	Each
Needle Holder, "V" Shape Self Righting, $\Phi 5 \times 330\text{mm}$	HK-101.024B	Each
Needle Holder, "V" Shape Self Righting, $\Phi 5 \times 330\text{mm}$	HK-101.024C	Each
Needle Holder, "V" Shape Self Righting, $\Phi 5 \times 330\text{mm}$	HK-101.024D	Each
Needle Holder, "V" Shape Straight Detachable, $\Phi 5 \times 330\text{mm}$	HK-101.024F	Each

### DISTRIBUTED BY

**Koala Medical Pty Ltd**

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: [sales@koalamedical.com.au](mailto:sales@koalamedical.com.au)

